



Ring lid&tube rack optional

**Specifications**

Model(integrated chamber)	TWB-11D	TWB-15D	TWB-22D	TWB-30D
Controller	PID			
Voltage	220V/50Hz 110V/60Hz			
Power(W)	800W	1000W	1200W	1500W
Tem-range(°C)	RT.+5~100°C			
Temperature resolution	0.1°C			
Temperature fluctuation	≤±0.5°C			
Division accuracy	≤±1°C			
Timer	999 minutes			
Recycle Design	Null			
Safety	Overcurrent protection; Over temperature alarm; Prevent dry burning function;			
Chamber Size(mm)	300*240*150mm	300*240*200mm	500*290*150mm	500*290*200mm
Product Size(mm)	340*280*240mm	340*280*290mm	550*350*240mm	550*350*290mm
Packing Size(mm)	390*330*290mm	390*330*340mm	600*400*290mm	600*400*340mm