## Description

380  $^{\circ}$ C general heating equipment. The device is widely used for heating in the school laboratory, industrial and mining , food processing, biochemistry, agriculture, petrochemicals engineering and pharmaceuticals industry for the liquid.

## Features

- 1.Once shaping shell is made of cold rolled steel sheet with surface Spraying treatment.
- 2. Nichrome heating wire, litter warming up time, long service life.
- 3. Easy operation, wide application.
- 4. Precise PID temperature control meets high accuracy experiment needs.
- 5."A" type are analog, "D" type are with digital display.
- 6."D" type sigle position mantles are with outer sensor, multi position mantles are with inner sensor, outer sensor optional.



## Heating Mantle (HM-A/D)

Туре	Analog							Digital						
Model	HM-100A	HM-250A	HM-500A	HM-1000A	HM-2000A	HM-3000A	HM-5000A	HM-100D	HM-250D	HM-500D	HM-1000D	HM-2000D	HM-3000D	HM-5000D
Capacity(ml)	100	250	500	1000	2000	3000	5000	100	250	500	1000	2000	3000	5000
Voltage(V)	200-240V/50-60Hz or 100-120V/50-60Hz							200-240V/50-60Hz or 100-120V/50-60Hz						
Max use temp(℃)	380							380						
Power(w)	100	150	250	350	450	600	800	100	150	250	350	450	600	800
Working time	continuous							continuous						
Exterior size(mm)	ф200*160			ф260х200	ф300х230	ф340х245	ф350х250	ф200x160		ф260х200	ф300х230	ф340х245	ф350х250	
Packing size(mm)	230×230×170			290x290x220	330x330x250	365x365x295	390x390x310	230x230x170			290x290x220		365x365x295	-
N/G.W.(kg)	2.5/2.8			4/6	5/7	7.5/8	8.5/9.2	2.5/2.8		4/6	5/7	7.5/8	8.5/9.2	